



# CATALOGUE

**MLUNGISI HEALTHCARE IS LEVEL ONE (1) CONTRIBUTOR TO B-BBEE.**

Tel: 011 025 1616 | Cell: +27 82 321 1191 | Email: [admin@mlungisihealthcare.com](mailto:admin@mlungisihealthcare.com)

Address: Suite A, Waterfall Crescent, Waterfall Office Park, Bekker Road, Midrand

[www.mlungisihealthcare.com](http://www.mlungisihealthcare.com)





## PLASMA CHEMICAL STRIP

### Application

Verify that the sterilization conditions are met during the sterilization process, such as temperature pressure, timing and other sterilization factors

### Features

1. The color block changes obviously before and after sterilization, which is easy to identify
2. Water-based ink, non-toxic, environmentally friendly
3. The indicator color printed firmly, won't melt and drop off during sterilization



## PLASMA TAPE

### Application

Sealing of flexible packaging materials such as non-woven fabrics, and visually display whether items have been sterilized.

### Features

1. The color contrast is clear and obvious before and after sterilization
2. Marking on the tape is OK
3. Easy to tear, without leaving any residual glue
4. The indicating color block is printed firmly fall off during sterilization.




## PLASMA TYVEK POUCH REEL

The packaging of metal, plastic and other medical materials sterilized by hydrogen peroxide.

### Features

1. Non-toxic water-based ink, color change clear and stable display, no melt.
2. Three explosion-proof edge technology can effectively prevent rupture
3. After hydrogen peroxide sterilization, the indicator blue color block changed pink.



## EO INDICATOR LABEL

### Application

1. ETO gas sterilization chemical indicator label sticker.
2. Red turn Blue when processed.
3. Verify that the items if got sterilized by ETO gas sterilizer.



## EO CHEMICAL STRIP

### Application

Verify that the sterilization conditions are met during the sterilization process, such as temperature pressure, timing and other sterilization factors

### Features

1. The color block changes obviously before and after sterilization, which is easy to identify
2. Water-based ink, non-toxic, environmentally friendly
3. The indicator color printed firmly, won't melt and drop off during sterilization



## EO INDICATOR TAPE

### Application

Sealing of flexible packaging materials such as non-woven fabrics, and visually display whether items have been sterilized.

### Features

1. The color contrast is clear and obvious before and after sterilization
2. Marking on the tape is OK
3. Easy to tear, without leaving any residual glue
4. The indicating color block is printed firmly fall off during sterilization.



## EO POUCH REEL

### Application

The packaging of metal, plastic and other medical materials sterilized by steam autoclave or EO.

### Features

1. Non-toxic water-based ink, color change clear and stable display, no melt.
2. Three explosion-proof edge technology can effectively prevent rupture
3. After pressure steam sterilization, the indicator blue color block changed black ; ethylene oxide sterilization, the red color block changed yellow.



## STEAM CHEMICAL STRIP

### Application

Verify that the sterilization conditions are met during the sterilization process, such as temperature pressure, timing and other sterilization factors

### Features

1. The color block changes obviously before and after sterilization, which is easy to identify
2. Water-based ink, non-toxic, environmentally friendly
3. The indicator color printed firmly, won't melt and drop off during sterilization



## STEAM POUCH REEL

### Application

The packaging of metal, plastic and other medical materials sterilized by steam autoclave or EO .

### Features

1. Non-toxic water-based ink, color change clear and stable display, no melt.
2. Three explosion-proof edge technology can effectively prevent rupture



#### Application

Sealing of flexible packaging materials such as non-woven fabrics, and visually display whether items have been sterilized.

#### Features

1. The color contrast is clear and obvious before and after sterilization
2. Marking on the tape is OK
3. Easy to tear, without leaving any residual glue
4. The indicating color block is printed firmly fall off during sterilization.

## STEAM CHEMICAL TAPE



## BD KIT

#### Application

Used in the measurement of the performance of the pre-vacuum pressure steam sterilizer

#### Features

1. Easy to mark: Use a marker pen to make relevant records on the test card. Such as: equipment number, operator, date, result, etc., easy to tracing .
2. The color change indication is clear, using high-temperature steam ink, the display is clear and stable, and will not melt



## LOW TEMPERATURE HYDROGEN PEROXIDE PLASMA STERILIZER



## LOW TEMPERATURE HYDROGEN PEROXIDE PLASMA STERILIZER




## LOW TEMPERATURE HYDROGEN PEROXIDE PLASMA STERILIZER



## **ETHYLENE OXIDE STERILIZER**



## **ETHYLENE OXIDE STERILIZER**



## **PRE-VACUUM PRESSURE STEAM STERILIZER**



## **PRE-VACUUM PRESSURE STEAM STERILIZER**



## **PULSATING VACUUM PRESSURE STEAM STERILIZER**



## **ENDOSCOPE CLEANING WORKSTATION**





## **AUTOMATIC WASHER DISINFECTOR**



## **AUTOMATIC WASHER DISINFECTOR**



## **PULSATING VACUUM PRESSURE STEAM STERILIZER DOUBLE DOOR**



## **WATER PURIFIER**



## **ULTRASONIC CLEANING MACHINE**



## **MEDICAL SEALING MACHINE**



## **EQUIPMENT TRANSFER VEHICLE**



## **STERILIZATION BASKET**



## **INSTRUMENT INSPECTION PACKING TABLE**

Tel: 011 025 1616 | Cell: +27 82 321 1191 | Email: [admin@mlungisihealthcare.com](mailto:admin@mlungisihealthcare.com)  
Address: Suite A, Waterfall Crescent, Waterfall Office Park, Bekker Road, Midrand

[www.mlungisihealthcare.com](http://www.mlungisihealthcare.com)



## DRESSING INSPECTION PACKING TABLE



**MLUNGISI HEALTHCARE IS LEVEL ONE (1) CONTRIBUTOR TO B-BBEE.**

Tel: 011 025 1616 | Cell: +27 82 321 1191 | Email: [admin@mlungisihealthcare.com](mailto:admin@mlungisihealthcare.com)  
Address: Suite A, Waterfall Crescent, Waterfall Office Park, Bekker Road, Midrand

[www.mlungisihealthcare.com](http://www.mlungisihealthcare.com)

